



Lab and ProLab Series

MEASUREMENT OF pH, ISE, CONDUCTIVITY AND
DISSOLVED OXYGEN - ACCURATE, RELIABLE AND SENSITIVE

SI Analytics

a xylem brand

The New Lab and ProLab Series

For all applications the right solution

SI Analytics is proud to present the new Lab and ProLab Series of pH, ISE, conductivity, and dissolved oxygen meters and accessories. The Lab and ProLab series of meters are ideal for any scientist looking to compliment their laboratory with reliable, robust and sensitive measurement readings.

For more information please visit our website www.si-analytics.com.

Lab 745, Lab 845 and Lab 945

User-friendly design for training and routine measurements.

Its intuitive operation and robust aluminum housing render the Lab x45 product series perfect for training and routine purposes.



Lab 845

Table of contents

Lab 745/845/945

Lab 855/865/955

IDS® Sensors

Lab 875 and 875P

ProLab 2500

ProLab 5000

Page 3

Page 6

Page 12

Page 16

Page 18

Page 22

Lab 845 Set



Lab 745 DO Meter, Lab 845 pH Meter,

Ordering information

| Type No. | Order No. | Description |
|--------------------|-----------|--|
| Lab 745 Set | 285206800 | Measuring ranges 0.01...20 mg/l DO. Set includes stand, power supply, and DO measuring cell Ox1113T |
| Lab 845 Set/BL19pH | 285206810 | Measuring pH, mV, ISE, temp., 3-point-cal., micropr., BNC connection. Set includes stand, power supply, BlueLine 19 pH, and DIN buffers in ampules (6 pieces) |
| Lab 845 Set/BL25pH | 285206820 | Measuring pH, mV, ISE, temp., 3-point-cal., micropr., BNC connection. Set includes stand, power supply, BlueLine 25 pH, and DIN buffers in ampules (6 pieces) |
| Lab 845 Set/BL29pH | 285206830 | Measuring pH, mV, ISE, temp., 3-point-cal., micropr., BNC connection. Set includes stand, power supply, BlueLine 29 pH, and DIN buffers in ampules (6 pieces) |
| Lab 945 Set/LF435T | 285206840 | Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temperature. Set includes stand, power supply, cond. cell LF435T, and cond. testing solution in ampules (6 pieces) |
| Lab 945 Set/LF513T | 285206850 | Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temperature. Set includes stand, power supply, cond. cell LF513T, and cond. testing solution in ampules (6 pieces) |
| Lab 945 Set/LF613T | 285206860 | Measuring ranges 0.000 µS/cm...500 mS/cm, salinity, total dissolved solids (TDS), temperature. Set includes stand, power supply, cond. cell LF613T, and cond. testing solution in ampules (6 pieces) |
| Z 611 | 285206380 | Connector, stand, and electrode holder for Lab 745/845/945 |
| Z 612 | 285206390 | Universal power supply for Lab 745/845/945 |
| Z 613 | 285206400 | USB cable with data transfer software for Lab 745/845/945 |
| Z 614 | 285206430 | Rubber pads for Lab 745/845/945 (4pcs) |
| Ox 1113T | 285206410 | Membrane covered amperometric sensor, plastic shaft, with temperature compensation, 1 m fixed cable with 8-pole plug, length 120 mm, 12 mm Ø, -5...+45 °C |
| LF 435T | 285206420 | 4 pole cell, plastic shaft, 1.5 m cable with 8 pole plug, sensor material graphite, cell constant 0.33 cm ⁻¹ , temp.-sensor NTC30kOhm, length 120 mm, 12 mm Ø, -5...+80 °C |
| Z 615 | 285206440 | Maintenance set for Ox1113T (3 x exchange heads, 10 x electrolytes) |
| Z 616 | 285206450 | Cable for connecting a RS232 printer to Lab 745/845/945 |

Technical data:

Lab 745 DO Meter

| | |
|--------------------------|--|
| Measuring range | 0 ... 200 %; 0 ... 20 mg/l; temperature: -10 ... 100 °C |
| Resolution | 1 %; 0.01 mg/l; 0.1 °C |
| Temperature compensation | Automatic with NTC30kOhm or fixed temperature |
| Accuracy | ±1 digit, ± 0.5 % of the measuring range, T [°C] ± 0.1 (5...50 °C) |
| Connectors | 8-pole sensor socket, 4-pole USB interface socket |
| Calibration | Direct input |
| | Temperature offset |
| | Single-point |
| | Automatic |
| Data storage | 4.000 Entries with date, time, value 1+2 and temperature |

Lab 945 Conductivity Meter

Lab 845 pH Meter

| | |
|--------------------------|--|
| Measuring range | pH: 0 ... 14; - 1,999 ... 1,999 mV; |
| | Temperature: -10 ... 100 °C |
| | ISE: 0 ... 30,000 ppm |
| Resolution | 0.01 pH; 1 mV; 0.1 °C |
| Accuracy | pH: ± 0.01 (± 2 pH around calibration point), U [mV] ± 0.3, T [°C] ± 0.1 (0..100 °C) |
| Temperature compensation | automatic with Pt1000 or fix temperature |
| Connectors | BNC, 2 x banana socket (4 mm), 4-pole USB Interface socket |
| Calibration | Direct input |
| | Temperature offset |
| | Three-point |
| | Automatic (DIN19266, TechDIN19267, Merck, Mettler) |
| Data storage | 4.000 Entries with date, time, value 1+2 and temperature |

Lab 945 Conductivity Meter

| | |
|--------------------------|--|
| Measuring range | 0 ... 200 µS/cm; 0 ... 2,000 µS/cm; 0 ... 20 mS/cm; 0 ... 500 mS/cm; |
| | Automatic range |
| | TDS: 0 ... 200 mg/l; 0 ... 2,000 mg/l; 0 ... 20 g/l; 0 ... 500 g/l |
| | Salinity: 0 ... 70 (after IOT) |
| | Temperature: -10 ... 100 °C |
| Resolution | 0.1 µS; 1 µS; 0.01 mS; 0.1 mS; 0.1 °C |
| Accuracy | ±1 digit, ± 0.5 % of the measuring range, T [°C] ± 0.1 (5...50 °C) |
| Temperature compensation | Automatic with NTC30kOhm or fixed temperature |
| Connectors | 8-pol sensor socket, 4-pole USB interface socket |
| Calibration | Direct input |
| | Temperature offset |
| | Single-Point |
| | Automatic |
| Data storage | 4.000 Entries with date, time, value 1+2 and temperature |

For all:

| | |
|---------------------|--|
| Display | Graphic LCD Display, 128 x 64 pixel, backlid |
| Interface | USB, isolated |
| Ambient temperature | -10 ... 55 °C |
| Housing protection | Aluminum IP40 |
| Dimensions | 145 x 185 x 55 mm (L x W x H) |
| Weight | Approximately 1 lb 9 oz (incl. power supply and stand) |
| EMC | Acc. EN 61326 class B |

Lab 855, Lab 865, and Lab 955

Precise. Reliable. Selective.

The Lab 855, Lab 865, and Lab 955 unite the most modern measuring technology available along with new functionality such as AutoRead and CMC (measuring range monitoring) which makes lab measurements even more reliable.

The newly designed, clearly structured keyboards are adapted to operators' logic with tactile feedback as well as large, easy-to-read displays which are used to support and enhance the interface between the meter and the user.



Precise measurements ...

... with Lab 855 and Lab 955



Reliable documentation ...

... with Lab 865



Precise measurements...

... with Lab 855 and Lab 955.



Modern meters for everybody who wants to simply measure accurately.

The Lab 855 for pH and Lab 955 for conductivity measurements are perfectly suited benchtop meters for measurements in laboratories in the chemical and pharmaceutical industries as well as in medical labs.


Precise measurements

Users achieve reproducible measured results due to the active automatic AutoRead function with independent detection of stable measuring values. An adjustable calibration timer assists in an increased improvement of the accuracy.

Easy to operate

The user-friendly keyboard with large, easy to read LCD display, deliver all relevant information at a glance.

| Type No. | Order No. | Description |
|-------------|-----------|--|
| Lab 855 Set | 285206700 | Simple, easy-to-use pH/mV benchtop meter (DIN) with universal power supply, stand and operating instructions, pH electrode BlueLine 14 pH, buffer solutions, 3 mol/l electrolyte solution. |
| Lab 955 Set | 285206760 | Simple, easy-to-use conductivity benchtop meter. Set includes conductivity measuring cell, device with universal power supply, stand, 4-pole graphite cell LF413T, and 0.01 mol/l KCl conductivity standard. |



- ▶ Reproducible measuring results with active AutoRead function
- ▶ Simple calibration with adjustable calibration timer
- ▶ Intuitive operation with clearly arranged keyboard

Benefits
Lab 855 / Lab 955

Also available as application-oriented sets with sensors, including power supply and stand.

Reliable documentation...

... with Lab 865



Precise measurements with documentation


The Lab 865 is perfect for pH measurements in quality assurance labs requiring the documentation for GLP. Built on the Lab 855 platform, the Lab 865 offers additional convenient functionality:

- **Data transfer via USB interface**
- **Manual or timer controlled data logging of measured values .**
- **Protocols include date, time, and ID and sensors serial numbers for GLP compliance.**
- **Data is transferred in *.csv format.**
- **Excel Add-in included for the formatted output of all data and calibration protocols.**

Also available as an application-oriented set with sensor, power supply and stand.

Lab 865

- USB interface for rapid data transfer



Data output in *.csv format



CMC function

- ▶ **Easy to use**
Graphic display with text menu for easy handling.
- ▶ **pH measurement on sight**
Optimize measuring results: With the new CMC function to monitor the congruency of measuring and calibration range for pH.

| Type No. | Order No. | Description |
|-------------|-----------|---|
| Lab 865 Set | 285206710 | Measuring parameters pH, mV, temp., 5-point-cal., micropr., Mini USB-B, data storage, DIN 19262 connect. Including stand, power supply, pH-temp. comb. electrode BlueLine 14 pH, calibr. solutions. |

Benefits
Lab 865

Lab 855, Lab 865, and Lab 955

Connectivity

Lab 855



Lab 865



Lab 955




Technical data

| Model | Lab 855 | Lab 865 | Lab 955 |
|--------------------------|--|--|---|
| Temperature compensation | Automatic/manual | Automatic/manual | Automatic, can be switched off |
| Calibration points | 1 to 3 | 1 to 5 pH, 2-7 ISE | 1 |
| Calibration records | 1 | 10 | 1 |
| Calibration timer | ■ | ■ | ■ |
| Memory entries | | 500/5000* | |
| Interface | | Mini USB-B | |
| GLP/AQS supporting | | ■ | |
| Display | LCD | Graphic b/w, backlit | LCD |
| Electrode connection | DIN | DIN | 8-pin |
| Additional | | CMC, input of sensor serial number | |
| Power supply | Battery or universal power supply | Battery or universal power supply | Battery or universal power supply |
| pH | - 2.0 ... 20.0 ± 0.1 pH - 2.00 ... 20.00 ± 0.01 pH - 2.000 ... 19.999 ± 0.005 pH | - 2.0 ... 20.0 ± 0.1 pH - 2.00 ... 20.00 ± 0.01 pH - 2.000 ... 19.999 ± 0.005 pH | |
| mV | ± 1200.0 ± 0.3 mV ± (2000 ± 1) mV | ± 1200.0 ± 0.3 mV ± (2500 ± 1) mV | |
| Temperature | - 5.0 ... 105.0 °C ± 0.1 °C | - 5.0 ... 105.0 °C ± 0.1 °C | |
| CMC | | ■ | |
| Conductivity | | | 0.00 ... 1000 mS/cm ± 0.5 % of meas. val. 0.000 ... 1.999 µS/cm, K = 0.01 cm ⁻¹ 0.00 ... 19,99 µS/cm, K = 0.1 cm ⁻¹ |
| Specific resistance | | | 0.00 ... 199.9 MΩcm |
| Cell constants fix | | | 0.01 cm ⁻¹ |
| with calibration | | | 0.450 ... 0.500 cm ⁻¹ 0.800 ... 0.880 cm ⁻¹ |
| adjustable | | | 0.090 ... 0.110 cm ⁻¹ 0.250 ... 2.500 cm ⁻¹ |
| Salinity | | | 0.0 ... 70.0 (nach IOT) |
| TDS | | | 1 ... 1999 mg/l |
| Temperature | | | -5.0 ... 105.0 °C ± 0.1 °C |
| T _{ref} | | | 20 °C/25 °C |
| Temperature compensation | | | none, nIF, 0.000 ... 3.000 %/K |

all measured values ± 1 decimal place

* manual/automatic

IDS®




New features

SI Analytics' IDS®: Intelligent, Digital Sensors technology for the standard parameters pH, conductivity and dissolved oxygen consists of two components, Digital sensors and matching field or benchtop meters. This new processing of the measured values no longer takes place in the device, exclusively in the sensor so that every sensor has its own data base when connected.

Built on the basic sensor of the BlueLine and ScienceLine series that have proven themselves tens of thousands of times over, the IDS® sensors have added precision and reliability and cover almost any application.

I intelligent:
IDS® sensors are intelligent. They log into the device automatically, submit their name, serial number, calibration status and history as well as all parameters.

D digital:
IDS® sensors transform the sensitive measuring signals in the sensor head into digital signals and transmit them to the output device without interference and errors.



S sensor:
IDS® sensors are based on proven and continuously developed sensors by SI Analytics. They cover almost any lab application, like pH, conductivity or dissolved oxygen measurements.



SI Analytics also offers Field meters with IDS®: **HandyLab 680**

IDS®-Sensors

Unique.

IDS® combines proven measuring technology with new advantages. Based on established electrochemical SI Analytics sensors, but equipped with state-of-the-art measuring electronics IDS® save the serial number and calibration data in the sensor, error-free and ready to use immediately. However, IDS® do not only store data, but they also process measuring signals directly and thus improve the data quality. This also allows a current evaluation of the sensor quality by means of the QSC (Quality Sensor Control) function.

IDS® combine proven technology with new advantages.

- High-quality, highly developed sensor technology combined with state-of-the-art measuring electronics.
- IDS® have saved the serial number and calibration history - error-free and therefore immediately ready for use.
- Current evaluation of the sensor quality for IDS® pH electrodes thanks to QSC (Quality Sensor Control).
- IDS® conductivity measurement: Two sensors to cover all applications.

SI Analytics
IDS®



IDS® types

| Type | Measuring function | Application range |
|-------------------------|--|--|
| A 157 IDS | pH micro cylinder + temp. | -5..100 °C pH 0-14 |
| A 161/162 IDS | pH sphere + temp. | -5..100 °C pH 0-14 |
| A 6880 IDS | pH spear + temp. | -5..100 °C pH 0-14 |
| A 7780 IDS | pH sphere + temp. | -5..80 °C pH 0-14 |
| BL 14 pH IDS | pH cone + temp. | -5..100 °C pH 0-14 |
| BL 21 pH IDS | pH spear | -5..80 °C pH 2-13 |
| BL 24/24-3 pH IDS | pH cylinder | -5..80 °C pH 0-14 |
| BL 27 pH IDS | pH flat | -5..50 °C pH 2-13 |
| BL 31 RX IDS | ORP platinum disk, 4 mm Ø + temp. | -5..100 °C |
| BL 32/32-3 RX IDS | ORP platinum pin, 1 mm Ø + temp. | -5..80 °C |
| FDO 1100/1100 3M IDS | Oxygen optical (photoluminescence) + temp. | 0..50 °C 0 ... 20 mg/l O ₂ |
| IL-Micro-pHT-IDS | pH micro cylinder + temp. | -5..100 °C pH 0..14 |
| IL-pHT-A120/ 170 MF-IDS | pH sphere + temp. | -5..100 °C pH 0..14 |
| IL-Sp-pHT-IDS | pH spear + temp. | -5..100 °C pH 0..14 |
| LF313T IDS | Conductivity stainless steel + temp. | -5..100 °C 0,01..200 µS/cm |
| LF 413T/413T 3M IDS | Conductivity graphite + temp. | -5..80 °C 1 µS/cm..2000 mS/cm |



- ▶ Higher accuracy than traditional analog sensors
- ▶ Perfect galvanized separation
- ▶ Resistant against environmental influences
- ▶ QSC takes the guess work out of the determining the health of your sensor
- ▶ Effortless capture and storage of your sensors latest calibration data
- ▶ Highest possible operator comfort and measuring precision

Benefits
IDS® Electrodes

Lab 875 and Lab 875P

Safe determination of pH, ORP, conductivity & optical DO...
... with the innovative Lab 875 and Lab 875P.

The Lab 875 with a digital measuring channel is optimal in the world of digital multi-parameter measurement using IDS®. The IDS® technology allows optimized measurements and efficient documentation in the simplest manner.

- One-channel multi-parameter meter for all IDS® sensors
- Digital sensor recognition
- Optionally installed printer: Lab 875P



Documentation as per GLP/AQS

- ▶ Automatic, digital capture of all IDS® sensor data for traceability of measured values.
- ▶ User administration capabilities for the safe allocation of user and measuring results.
- ▶ Transmission of all data in *.csv format via USB interface to PC; if desired, formatted transfer to Excel (MultiLab® Importer, included in the delivery and as a download).
- ▶ Output directly into an optional integrated printer.


Benefits
Lab 875 and Lab 875P

ProLab 2500

If there is the need for a little more...

The ProLab 2500 is a digital high-performance meter with three channels. With its color graphic display protected by glass, high-quality zinc die-casting base as well as the anti-bacterial keyboard cover, it meets even the highest demands.

- Three universal measuring channels
- Digital sensor recognition
- Antibacterial keyboard
- Any configuration of pH, ORP, conductivity and Optical DO



ProLab 2500



Flexible performance



- Measures pH, ORP, ISE, dissolved oxygen, and conductivity.
- Any combination of the same or different parameters.
- Color backlit graphic display with CMC, QSC and channel display.
- An adapter for conventional pH electrodes (DIN plug) included in the delivery.
- Memory with 10,000 possible entries.

Measuring certainty

- The digital signal transfer eliminates interference, safely allocates calibration data, automatically transmits sensor data.
- The Quality Sensor Control (QSC) icon provides precise information about the actual state of the electrode and therefore increases the operational safety.

Documentation

- Automatic, digital capture of all sensor data for retraceability of measured values.
- Available user administration for the safe allocation of user and measuring results or sample and measuring results.
- Data output to PC, USB memory stick or selected printers.



Digital sensor recognition



Lab 875 (P) and ProLab 2500

Technical data

| Model | Lab 875 (P) | ProLab 2500 |
|------------------------------------|---|--|
| Parameter | pH, mV, saturation, concentration, partial pressure, conductivity, spec. resistance, salinity, TDS, temperature | |
| IDS® - Intelligent Digital Sensors | ■ | ■ |
| Universal measuring channels | 1 | 3 |
| Analog pH/ORP sensors | Z600 IDS® analog adapter (optional) | Z600 IDS® analog adapter (included in delivery) |
| Temperature compensation | all except for ORP | |
| Calibration points | 1-5 | 1-5 |
| pH | - | 2-7 (requires adapter) |
| ISE | 1 | 1 |
| Dissolved oxygen | 1 | 1 |
| Conductivity | 1 | 1 |
| Calibration records | Max. 10 | Max. 10 |
| Calibration timer | 1 - 999 days | 1 - 999 days |
| Memory entries | manual: 500 data sets automatic: 4,500 data sets | manual: 500 data sets automatic: 10,000 data sets |
| Logger | ■ | ■ |
| Interface | Mini USB-B | USB-A, Mini USB-B |
| GLP/AQS supporting | ■ | ■ |
| Display | Graphic, BW | Color graphic |
| Printer option | Yes: Lab 875P | external |
| Additional | CMC, QSC | antibacterial keypad, QSC, CMC, replaceable firmware |
| Power supply | Universal power supply, battery (4 x 1.5 V AA Type) | Universal power supply |



ProLab 2500

Lab 875P



Ordering information

| Type No. | Order No. | Description |
|----------------------------|-----------|---|
| Lab 875 | 285206320 | One channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand and power supply. |
| Lab 875 pH Set | 285206720 | One channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand, power supply, BlueLine 14 pH IDS® and buffer solutions. |
| Lab 875 Cond Set | 285206730 | One channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand, power supply, LF 413T IDS® and cond. testing solutions. |
| Lab 875P | 285206330 | One channel instrument with integrated printer for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand and power supply. |
| Lab 875P pH Set | 285206740 | One channel instrument with integrated printer for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand, power supply, BlueLine 14 pH IDS® and buffer solutions. |
| Lab 875P Cond Set | 285206750 | One channel instrument with integrated printer for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., Mini USB-B interface, Data storage. Including stand, power supply, LF 413T IDS® and cond. testing solutions. |
| ProLab 2500 | 285206350 | Three channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., USB-A and Mini USB-B interface, Data storage. Including stand, power supply and IDS® analog adapter. |
| ProLab 2500 pH Set | 285206770 | Three channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., USB-A and Mini USB-B interface, Data storage. Including stand, power supply, A 162 IDS®, buffer solutions and IDS® analog adapter. |
| ProLab 2500 pH/Cond Set | 285206780 | Three channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., USB-A and Mini USB-B interface, Data storage. Including stand, power supply, A 162 IDS®, LF 413T IDS®, testing solutions and IDS® analog adapter. |
| ProLab 2500 pH/Cond/Ox Set | 285206790 | Three channel instrument for IDS® sensors in order to measure pH, mV, Cond., DO and Temp., USB-A and Mini USB-B interface, Data storage. Including stand, power supply, A 162 IDS®, LF 413T IDS®, FDO 1100 IDS®, testing solutions and IDS® analog adapter. |
| Z 600 | 285206360 | Adapter IDS socket / DIN plug for connecting analog DIN electrodes to an IDS® socket. |
| Z 610 | 285206370 | Printer paper, document quality, one roll for Lab 875P. |
| Z 850 | 285204889 | Universal power supply unit, 230 and 120 V for all Lab and ProLab meters. |
| Z 865 | 285201520 | Stand set S4D, including arm and electrode holder for docking to the meters of the Lab- and ProLab family as well as for autonomous usage. |
| Z 866 | 285204940 | Flexible electrode arm for fixed attachment to Lab & ProLab meter family. |
| Z 875 | 285201540 | USB cable for Lab and ProLab meters. |



ProLab 5000

Complete system for measurement, control, and regulation of pH, conductivity, and dissolved oxygen in laboratories and technical centers.

- ▶ Measurement of pH/ISE/mV, conductivity, and dissolved oxygen
- ▶ Up to 4 measurement modules (inputs) in a variety of configurations
- ▶ 5.7" RGB TFT display with touch control
- ▶ PC software with extensive operating functions
- ▶ Coupling of auto-sampler and burettes for dosing and automated measurements
- ▶ Special electrodes for pH
- ▶ Current outputs for each parameter
- ▶ Additional modules for current output available
- ▶ Timer function
- ▶ Alarm/threshold function
- ▶ 2 PID regulators
- ▶ Virtual channels to calculate different parameters from the measured value

- ▶ Data storage and recording
- ▶ Data transfer with RS232/USB or ethernet
- ▶ Logbook (i.e. documentation of setting changes)
- ▶ Password protection

ProLab 5000




Components of the ProLab 5000

The ProLab 5000 multiparameter measuring instrument is based on a modular system structure that links the measurement modules to the central unit and to each other via a bus connection. It is the complete solution for measurement, control, and regulation in laboratory and technical centers. Choose up to four measurement modules can determine any combination of pH, conductivity, dissolved oxygen, ORP/redox potential, temperature and other parameters.

The main components of the multiparameter laboratory measurement system ProLab 5000 are:

- ▶ ProLab 5000 with power supply, touch screen display, internal modules, data logger and logbook, various digital interfaces such as RS232, USB and ethernet
- ▶ PC visualization program
- ▶ Measurement modules
- ▶ Sensors
- ▶ Optional samplers, depending on type, up to 72 samples (i.e. TW Alpha plus)
- ▶ Optional dosing system (TITRONIC® 500)
- ▶ Optional stirrer

Setting options



Multiparameter laboratory measuring instrument ProLab 5000

Color graphical touchscreen display

The distinguishing outward feature of the ProLab 5000 is the color touchscreen graphical display, which significantly simplifies operation, configuration, and calibration. Up to four measured values with their units, temperature values, a possible second value as well as an individual name can be presented simultaneously on a display page.

Four integrated threshold relays

Four integrated threshold relays for simple control, regulation, or alarm tasks are standard. Two PID-regulators that work independently of each other are available for complex regulation tasks. The regulators can be assigned to any parameter and function as analog, pulse width or pulse frequency regulators using the analog current outputs and the relay outputs of the measuring instrument.

ProLab 5000 Pilot visualization program

All values can be cyclically output and displayed graphically via the USB interface or ethernet connected to the PC visualization program ProLab 5000 Pilot. The visualization program ProLab 5000 Pilot simultaneously contains all functions for reading the data logger. The data logger can be used via the USB interface and the data can be transferred directly to a USB without the program ProLab 5000 Pilot.

Data recorder

The data recorder graphically displays the measured value curve of the parameters of each measurement module over the past 48 hours on the touch screen. This representation gives the operator a quick overview of the measurement progress, indicating the average, maximum, and minimum readings.



ProLab 5000

Ordering information

| Type No. | Order No. | Product hierarchy No. | Description |
|-----------------------|-----------|-----------------------|--|
| PL5000 0D 1pH 0LF 0OX | 285206010 | 07 | Multiparameter measuring instrument ProLab 5000 with one pH/mV/ISE module |
| PL5000 0D 1pH 1LF 0OX | 285206020 | 07 | Multiparameter measuring instrument ProLab 5000 with one each pH/mV/ISE and conductivity module |
| PL5000 0D 1pH 1LF 1OX | 285206030 | 07 | Multiparameter measuring instrument ProLab 5000 with one each pH/mV/ISE, conductivity and dissolved oxygen module |
| PL5000 0D 1pH 2LF 1OX | 285206040 | 07 | Multiparameter measuring instrument ProLab 5000 with one each pH/mV/ISE and dissolved oxygen as well as two conductivity modules |
| PL5000 0D 2pH 0LF 0OX | 285206050 | 07 | Multiparameter measuring instrument ProLab 5000 with two pH/mV/ISE modules |
| PL5000 0D 2pH 1LF 0OX | 285206060 | 07 | Multiparameter measuring instrument ProLab 5000 with two pH/mV/ISE and one conductivity modules |
| PL5000 0D 2pH 2LF 0OX | 285206070 | 07 | Multiparameter measuring instrument ProLab 5000 with each two pH/mV/ISE and conductivity modules |
| PL5000 0D 2pH 1LF 1OX | 285206080 | 07 | Multiparameter measuring instrument ProLab 5000 with one each conductivity, dissolved oxygen and two pH/mV/ISE modules |
| PL5000 0D 3pH 0LF 0OX | 285206090 | 07 | Multiparameter measuring instrument ProLab 5000 with three pH/mV/ISE modules |
| PL5000 0D 3pH 1LF 0OX | 285206100 | 07 | Multiparameter measuring instrument ProLab 5000 with three pH/mV/ISE and one conductivity modules |
| PL5000 0D 3pH 0LF 1OX | 285206110 | 07 | Multiparameter measuring instrument ProLab 5000 with three pH/mV/ISE and one dissolved oxygen modules |
| PL5000 0D 4pH 0LF 0OX | 285206120 | 07 | Multiparameter measuring instrument ProLab 5000 with four pH/mV/ISE modules |
| PL5000 1D 1pH 0LF 0OX | 285206130 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter and one pH/mV/ISE module |
| PL5000 1D 1pH 1LF 0OX | 285206140 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter as well as one each pH/mV/ISE and conductivity module |
| PL5000 1D 1pH 1LF 1OX | 285206150 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter as well as one each pH/mV/ISE, conductivity and dissolved oxygen module |
| PL5000 1D 1pH 2LF 1OX | 285206160 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter as well as one each pH/mV/ISE, dissolved oxygen and two conductivity modules |
| PL5000 1D 2pH 0LF 0OX | 285206170 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog-converter and two pH/mV/ISE-modules |
| PL5000 1D 2pH 1LF 0OX | 285206180 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog-converter, two pH/mV/ISE- and one conductivity-modules |
| PL5000 1D 2pH 2LF 0OX | 285206190 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog-converter, each two pH/mV/ISE- and conductivity-modules |
| PL5000 1D 2pH 1LF 1OX | 285206200 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter, one each conductivity and dissolved oxygen as well as two pH/mV/ISE modules |
| PL5000 1D 3pH 0LF 0OX | 285206210 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter and three pH/mV/ISE modules |
| PL5000 1D 3pH 1LF 0OX | 285206220 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter, three pH/mV/ISE and one conductivity modules |
| PL5000 1D 3pH 0LF 1OX | 285206230 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter, three pH/mV/ISE and one dissolved oxygen modules |
| PL5000 1D 4pH 0LF 0OX | 285206240 | 07 | Multiparameter measuring instrument ProLab 5000 with digital-analog converter and four pH/mV/ISE modules |
| Z570 | 285206250 | 04 | USB cable for ProLab 5000 |
| Z573 | 285206260 | 04 | Cable for connecting the ProLab 5000 to sample changer |
| Z575 | 285206270 | 04 | Wall mount for ProLab 5000 |


| pH/mV/ISE | Number of modules | | |
|-----------|-------------------|----|--------------------------|
| | Cond. | DO | Digital-analog converter |
| 1 | 0 | 0 | 0 |
| 1 | 1 | 0 | 0 |
| 1 | 1 | 1 | 0 |
| 1 | 2 | 1 | 0 |
| 2 | 0 | 0 | 0 |
| 2 | 1 | 0 | 0 |
| 2 | 2 | 0 | 0 |
| 2 | 1 | 1 | 0 |
| 3 | 0 | 0 | 0 |
| 3 | 1 | 0 | 0 |
| 3 | 0 | 1 | 0 |
| 4 | 0 | 0 | 0 |
| 1 | 0 | 0 | 0 |
| 1 | 1 | 0 | 0 |
| 1 | 1 | 1 | 1 |
| 1 | 2 | 1 | 1 |
| 2 | 0 | 0 | 1 |
| 2 | 1 | 0 | 1 |
| 2 | 2 | 0 | 1 |
| 2 | 1 | 1 | 1 |
| 3 | 0 | 0 | 1 |
| 3 | 1 | 0 | 1 |
| 3 | 0 | 1 | 1 |
| 4 | 0 | 0 | 1 |



Technical specifications

| | |
|-------------------------------|---|
| Auxiliary energy | 24 V DC desktop power supply, voltage disconnect via the unit switch |
| Ambient temperature | 0 ... + 40 °C |
| Display | Touch screen graphic display 320 x 240 pixel, 256-color, back lit |
| Menu languages | German, English |
| Data transmission | Serial interface RS-232, ethernet port, USB port for PC Connection |
| Control outputs | 4 potential-free relay outputs; resistive load $I \leq 1$ A, $U \leq 24$ V DC for threshold or alarm function; including a relay with timer function (wash contact; adjustable time interval 1 ... 9,999 hours) |
| Data storage | Integrated data logger for approx. 100,000 values including date and time, 48 hour data recorder |
| Logbook | Approx. 200 activities including date and time |
| Housing | Aluminum housing IP 40/DIN EN 60529; Dimensions see dimensional drawings, wall-mount kit available |
| Connections | BNC, banana, 8-pin Din, BK, USB, ethernet |
| Electromagnetic compatibility | 89/336/EEC, EN 61326 Class B |
| Measurement modules | Four internal measurement modules; in any combination; inputs galvanically separated; calibration data storage; sensor monitoring via adjustable threshold bands; manual and automatic temperature compensation |
| Regulator module | Optional regulator module PL5000DAC: Standard signal module 4 x 0(4) ... 20 mA |
| Safety | Protection class III, EC Directive 73/23 EN 61010-1: 2001 |
| GLP | GLP functions (data recording) |

ProLab 5000 - dimensions:



Multifunctional connections:



Back panel connections

Lateral connections



Measurement modules:

| Internal modules | Main parameter measuring range/resolution | Secondary parameters/resolution | Temperature measurement measuring range/resolution | Electrodes/sensors |
|---------------------|--|--|--|---|
| PL5000 pH; ORP, ISE | pH value pH 0 ... 14 / 0.001 pH Accuracy: 0.005 pH | Chain voltage 0.1 mV | -10...130 °C / 0.1 °C | pH single rod measuring cell, separated measuring chain, Pt 1000 temperature sensor |
| | Redox potential -2000... +2000 mV > 0.1 mV | Redox voltage relative to the standard hydrogen electrode 0.1 mV | -10...130 °C / 0.1 °C | Redox single rod measuring cell, separate measuring chain, Pt 1000 temperature sensor |
| | Ion concentration corr. sensor specification (ISE) and calibration | Chain voltage 0.1 mV | -10...130 °C / 0.1 °C | Ion-selective electrode (ISE), separate measuring chain, Pt 1000 temperature sensor |
| PL5000LF | Conductivity 0...200 µS/cm 0.1 µS/cm (4-pol) 0...2 mS/cm 1 µS/cm 0...20 mS/cm 0.01 mS/cm 0...500 mS/cm 0.1 mS/cm automatic switching | Salinity 2 ... 42 g/kg | -10...130 °C / 0.1 °C | 4-electrode conductivity measurement cell, NTC30kOhm temperature sensor |
| PL5000OX | O ₂ saturation 0...120% / 0.1% | O ₂ concentration 0...20 mg/l / 0.1 | -10...130 °C / 0.1 °C | Membrane-covered amperometric O ₂ sensor, NTC30kOhm temperature sensor |

We are Xylem Analytics

Xylem consists of three business sectors - Water Solutions, Applied Water Systems and Analytics. The following brands make up Xylem Analytics and act like SI Analytics in the chemical, pharmaceutical, biotechnological, food and plastics industries.

Bellingham & Stanley

For almost a century, Bellingham + Stanley has been a well-known leader in the field of three product groups.

Our products include:

- Refractometers
- Polarimeters
- Certified Reference Materials

www.bellinghamandstanley.com



a xylem brand

SI Analytics

The brand SI Analytics offers nearly 80 years of experience in glass technology and the development of analysis equipment. We continue to develop and manufacture our products with the highest standards to fuel the requirements of innovation and quality our customers have come to expect. We continue to value tradition and manufacture in the footsteps of time-honored Mainz glass manufacturers.

Our products include:

- Electrodes
- Meters
- Titrators
- Capillary viscosimeters

www.si-analytics.com

SI Analytics

a xylem brand

ebro

Our brand ebro has been servicing the scientific world with innovative temperature measurement solutions for over forty years. Customer feedback has continued to play an important role within the business model. To ebro, customer care not only means supporting existing product and software, it also means being able to provide custom solutions enabling our customers to ensure their success.

Our products include:

- Precision thermometers
- Food Safety test kits
- Frying oil monitors
- Humidity, vacuum & temperature dataloggers
- Portable digital refractometers

www.ebro.com



OI Analytical

Since 1963 OI Analytical has been providing innovative products used for chemical analysis and is a key supplier of sample preparation and turn-key analytical solutions for testing food products and water for chemical contaminants.

Our products include:

- Total Organic Carbon (TOC) & cyanide Analyzers
- Automated Chemistry Analyzers
- GC Continuous Monitoring Systems
- Refrigerant Monitors
- Process/On-line Instruments
- GC Detectors and Systems

www.oico.com



Furthermore, Xylem Analytics comprises:



a xylem brand



a xylem brand



a xylem brand

Xylem | 'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water technology company.

We're a global team unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to www.xyleminc.com

SI Analytics
a **xylem** brand

**Xylem Analytics Germany
Sales GmbH & Co. KG**

SI Analytics
Hattenbergstr. 10
55122 Mainz
Germany

Phone: +49.6131.66.5111
Fax: +49.6131.66.5001
E-Mail: si-analytics@xyleminc.com
Internet: www.si-analytics.com

In North America:

P.O. Box 9010
151 Graham Road
College Station, Texas 77842-9010, USA

Toll-free: 866-691-7954
Local: 979 690 5563
E-Mail: information.request@xyleminc.com
Internet: www.si-analytics.com

presented by

SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.

© 2016 Xylem, Inc. 980 093US Version 09/2016