

Application

**Determination of betaine, alkyl-
betaine, cocamidopropyl-
betaine with perchloric acid**



Application

Appliances

- Titrator: TL 7000

Electrodes

- Electrodes: N 6480 Eth

Reagents

- Titrant: Perchloric acid (HClO_4) = 0.1 mol/l in 1,4-Dioxane
(in a 500 ml volumetric flask with 300 ml Dioxan add slowly 5,6 ml Perchloric acid 60 % p.a. and fill up with Dioxane to 500 ml. The solution is standardized with Kaliumhydrogenphthalate in 60 ml Acetic acid similar to HClO_4 in glacial acetic acid)
- Methanol p.a.
- Methyl-glycol p.a. (Ethylene glycol monomethyl ether)
- Basification solution (4 g Sodium hydroxide p.a. and 8 g Sodium acetate p.a. dissolved and fill up to 100 ml dest. water)

Use

The method describes the non-aqueous titration of Betaine in methanol/methyl-glycol with perchloric acid 0.1 mol/l in Dioxane. The sample is dissolved in methanol and basified with NaOH/NaAcetate solution to deprotonating all acids compounds contain the sample. After a reaction time the methyl glycol is added and the sample is titrated with HClO_4 in Dioxane. The complete curve give three EQs:

EQ1 is the NaOH excess,

EQ2 is the Na-Acetate and acids or amines contains the sample,

EQ 3 is the Betaine. The Betain concentration would them calculated EQ3-EQ2.

The first EQ is for the Betaine determination not important. The method includes a pretitration of 6 ml to negate the first EQ (NaOH)! This simplifies the evaluation.

The rest of the titration curve shows two EQs, the Betaine is calculated as the difference EQ2-EQ1.

Description

For Samples with 30 % (BASF I):

Into a 150 ml baker weigh 0,6000 g Betaine sample, dissolve in 20 ml methanol, added 0,5 ml basification solution and stir for 5 to 8 minutes. Then add 20 ml methanol and 60 ml methyl-glycol. The sample is titrated with HClO_4 0.1 mol/l in Dioxane. The Method includes a pretitration of 6 ml with 30 sec waiting time, it's important to overtitrate the first EQ (NaOH) for a correct calculation. The electrode is rinsed first with dest. water then with alcohol before a next and first titration.

For Samples with 1 % (BASF II):

Into a 150 ml baker weigh 5,0000 g Betaine sample, dissolve in 40 ml methanol, add 0,5 ml basification solution and stir for 5 to 8 minutes. Then add 60 ml Methyl glycol. Titrate the Sample with HClO_4 0.1 mol/l in Dioxane. The Method includes a pretitration of 6 ml with 30 sec waiting time, it's important to overtitrate the first EQ (NaOH) for a correct calculation. The electrode is rinsed first with dest. water then with alcohol before a next and first titration.

Application


Examples / Results

BASF (I)			
Weight in g	EQ1 in ml	EQ 2 in ml	Betaine in %
0,6047	10,532	15,590	29,72
0,6201	10,518	15,713	29,77
0,6069	10,542	15,615	29,70
		Mean.	29,73
		RSD in %	0,12


BASF (II)			
Weight in g	EQ1 in ml	EQ 2 in ml	Betaine in %
5,0301	11,294	13,1	1,28
5,0961	11,284	13,107	1,27
5,0661	11,381	13,131	1,23
		Mean.	1,26
		RSD in %	2,10

Application

BASF I (complete curve):



BASF I (method with 6 ml pretitration):



Method data

Method name:	Betaine with HClO ₄	Titration duration:	15 m 29 s
End date:	23.06.16	End time:	17:43:54

Titration data


Sample ID:	BASF 1	Weight:	0.62010 g
Start mV (A):	-175.5 mV	End mV (A):	499.8 mV

EQ1:	10.518 ml / 140.5 mV	EQ1:	10.518 ml
------	----------------------	------	-----------


EQ2:	15.713 ml / 388.0 mV	Betaine:	29.77 %
------	----------------------	----------	---------

Application

BASF II (complete curve):



BASF II (method with 6 ml pretitration):



Method data

Method name:	Betaine with HClO ₄	Titration duration:	10 m 46 s
End date:	24.06.16	End time:	16:33:42

Titration data

Sample ID:	BASF 2	Weight:	5.09610 g
Start mV (A):	-183.5 mV	End mV (A):	434.2 mV

EQ1:	11.284 ml / 188.9 mV	EQ1:	11.284 ml
------	----------------------	------	-----------

EQ2:	13.107 ml / 323.3 mV	Betaine:	1.27 %
------	----------------------	----------	--------

Application

Hints

SI Analytics GmbH
Hattenbergstr. 10
55122 Mainz
Germany
Homepage:

Phone: +49 (0) 6131 / 66 – 5062
+49 (0) 6131 / 66 – 5118
Fax: +49 (0) 6131 / 66 – 5001
E-Mail: titration@si-analytics.com
www.si-analytics.com